


PSS

PST

 **Precise**
Implant Systems

PHI

PSI



Catalog 2014



Precise



Table of contents:

Company Profile _____	4-5
Implant System _____	6
Advanced Implants _____	8 - 13
Dynamic Implants _____	14 - 20
Implant Surface _____	20 - 21
Implant Packaging _____	22 - 23
Surgical Kit _____	24 - 26
Surgical Drilling Sequences _____	27
Prosthetics _____	28
Healing Caps _____	30
Impression Coping _____	32
Peek Abutments _____	34
Straight Titanium Abutments _____	36
Angulated Titanium Abutments _____	38
Castable Abutments _____	40
Multi Unit Abutments _____	42
Ball Abutments _____	44



■ Company profile

We are the same but much better!

Precise company believes in making dental implants by redefining precision in the process of manufacturing, marketing and service. We combine this new level of precision with simplicity and quality control together with state of the art scientific knowledge, to achieve vast improvements over existing implants.

The partners in Precise have more than 30 years of experience as leading developers and marketers of implants and prosthetic parts, with proven sales of more than a million implants by other brands.

Precision

In Precise, we focus on making precision implantology. The products that we manufacture are of the highest quality, and they provide easy-to-use solutions. We strive for total customer satisfaction by producing the ultimate kit that addresses any clinical need.

Science

Precise is based on the principle of science by nature. This principle is carried out to its practical applications, which reflect our solution-driven technological innovations.

High Standards

Precise is committed to the highest level of accuracy and quality control to produce the best possible implants. The company's products comply with all the requirements of medical devices and are 1183 CE approved .

Product Range

Precise manufactures four types of dental implants and a variety of prosthetic parts, as well as extremely accurate, internal-hex, "one platform for all diameters" components, which enable the dentist to maintain a minimal inventory of prosthetic parts.

Our products are designed for maximum bone-to-implant contact (BIC), maximum initial stability, safe and easy insertion and high aesthetics.



Precise Implants

The geometry and features of our family of dental implants, along with our meticulous manufacturing process, result in implants that ensure maximum precision, smooth insertion and a more comfortable, natural healing process for the patient.

The implant family uses classic implant design and surface treatments, and an intermediate micro-thread, to achieve larger bone-to-implant surface area contact.

Precise Implants are made of biocompatible, extremely strong Ti-6Al-4V ELI medical-grade titanium alloy. Features such as flutes cut at the apex and a conical shape render the implant self drilling and self tapping. These make the implant insertion process as streamlined as possible, resulting in firmer contact with the bone.



PSI Precise Spiral Implant



PST Precise Spiral Top



PSS Precise Spiral Standard





PHI Precise Helical Implant

Precise Nano-Mesh Surface

Our implants are treated with a state-of-the-art process and high purity surface that promotes bone growth and accelerating osseointegration.

A complex process of blasting biocompatible particles followed by dual thermal acid etching.



Precise Nano-Mesh characteristics are:

- Hybrid macro and micro surface structure
- High oxide layer
- Excellent bio-compatibility surface
- Faster Bone-To-Implant-Contact (BIC)
- High wettability characteristics



■ Surgical instruments

Our surgical kit is comprised of required basic surgical instrumentation for implantation procedures



P-SKP 94069
Surgery Kit Plastic

P-SKS 94068
Surgery Kit Small SSL

■ Surgical Kit

The Precise surgical kit consists of basic surgical instrumentation for carrying out all implantation procedures.

Autoclaveable kit, user-friendly and simple to use, with optimal and convenient tool arrangement.

- Ergonomically box design
- Compact and light to carry
- Shock-resistant plastic compounds certified for over 1000 autoclave sterilization cycles
- Carefully arranged and sectioned tray to accommodate instruments in a separate compartment
- Individual silicon fixtures that holds instrument and prevent movement during transport, even if the kit is turned upside down
- Ware resistant and easy to clean Laser-etched markings
- Dimensions: 18.9 x 13.9 x 6.1 cm

Materials:

- Box: Radel R polyphenylsulfone resin
- Tool Holders: Autoclavable medical-grade silicon
- Bath: Stainless steel

Four different working positions with non-skid tray

■ Precise Prosthetics

Precise's diverse prosthetic array provides implant professionals in the dental industry with comprehensive solutions for all the common dental restorative procedures.

The meticulously designed prosthetic system features straight, angled and casting abutments, as well as ball attachments in various lengths, permitting versatile restoration with narrow implants. Our unique single restoration platform enables the use of any abutment of any implant diameter.



■ surgical drills



Precise surgical drills are available in variety of diameters and are made of surgical stainless steel. The drills are color coded easy diameter identification and Laser-marks on depth grooves for easy reference during surgery.



■ Conical Implants

The conical spiral implants are the heart of the implant world. The implants are ideal for soft bone of types III and IV, and even for softer types. Moreover, in special cases they may be used also in hard bone* of type I and II.



■ PSI
Precise Spiral Implant

With its science, design and quality control, this implant introduces a whole new level of stability, comfort and bone completion. This is why we can proudly assert that Precise applies the principle of science by nature, carrying out its technological innovation to practical solutions.



■ PST
Precise Spiral Top

Based on PSI's characteristics and advantages, this implant is at the top of the dental implant world today in terms of aesthetics.



■ PSS
Precise Spiral Standard

An improved iteration of the standard spiral conical implant with superior finishing, which greatly outperforms competing products.

* This does not constitute a recommendation for the dentist on how to use these implants.



■ Cylindrical Implants

■ PHI
Precise Helical Implant

Our classical, mainstream, implant is an indispensable element of every implant system available in the market. The PHI promotes bone condensing and creates overall stability, it is ideal for hard bone but suitable for common cases too, and it enables small changes in the implantation direction.



■ PSI ■ PSS ■ PST ■ PHI

The advantages of using Precise's conical implants:


- These implants enable penetration into locations with thicker diameter and in soft tissue and softer bones (bone types III and IV)
- The narrow conical implant is suitable for cases of limited bone width and it prevents damage to the anatomical structure
- The implants enable reasonable alteration of the implantation direction during the implant procedure
- These implants are ideal for immediate implantation as well as for immediate loading
- They can be indicated for any condition, bone type and surgical protocol
- They promote bone condensing and create enhanced primary stability
- The threads in the upper (coronal) part are thicker square, in the middle part are thinner square and in the apical part, where the borders are straight, the threads are sharper and V-shaped

To summarize, the advantages mentioned herein constitute an implant that enables easy insertion even to narrow locations with high primary stability, prevents damage to the anatomical structure and is self tapping.

The clinical advantages:

Our implants allow for better bone condensing and offer high primary stability and minimum bone loss, reducing trauma. They can be implanted in areas with small diameter, preventing damage to adjacent teeth. They can be precisely controlled during the implantation through the insertion pathway, meaning that their angle can be changed during the procedure to reduce the perforation in the lingual and buccal cortex.

The implants' enhanced stability is achieved even with small bone mass and in soft bone areas – necessary in many clinical cases. This is the result of the special implant body, the thread design and the apical part.





PSI

■ PSI
Precise Spiral Implant

With its science, design and quality control, this implant introduces a whole new level of stability, comfort and bone completion. This is why we can proudly assert that Precise applies the principle of science by nature, carrying out its technological innovation to practical solutions.



SECTION A-A
SCALE 7:1

NOTE:

① - Critical Dimensions

NORIS Medical

IMPLANT HO



■ PSI

Precise Spiral Implant




The next generation of the conical implants!

The PSI implant has a complex, new differential thread form that becomes more pronounced the more it progresses to the top of the implant. This structure features two tunnels that provide better fixation to the bone by at least 30% compared to implants with standard thread forms. The differential fixation channels enable much easier insertion toward the top, to achieve the excellent stabilization familiar to anyone working with conical implants.



Internal Hex connection

Is extremely precise, being the platform for all diameters that enables a simple restoration process.



Upper Part

Surface is rough towards the top of the implant, reduces bone resorption and reduces crestal stress for better load distribution.

2 Micro Rings



Implant Structure

Is conical, to achieve better bone condensing, higher initial stability and easy insertion.

Implant Threads

Are double threaded (2x2.4mm) for efficient implantation. The threads in the upper part are thicker square, in the middle part thinner square and in the lower part V-shaped to enable fast insertion, bone condensing and self tapping.




Lower part

Has narrow, deep and sharp apical starting threads that enable dual entrance, self tapping and self drilling, as well as easy bone penetration to prevent damage to the anatomical structure.






■ PSI



D1



PSI 3.3

Dia.	3.3mm	Dia.	3.3mm	Dia.	3.3mm	Dia.	3.3mm	Dia.	3.3mm
L	8.0mm	L	10.0mm	L	11.5mm	L	13.0mm	L	16.0mm
Ref.	91001	Ref.	91002	Ref.	91003	Ref.	91004	Ref.	91005
D1	3.70	3.70	3.70	3.70	3.70	3.70	3.70	3.70	3.70
D2	2.95	2.95	2.95	2.95	2.95	2.95	2.95	2.95	2.95
D3	3.50	3.50	3.50	3.50	3.50	3.50	3.50	3.50	3.50

PSI 3.75

Dia.	3.75mm	Dia.	3.75mm	Dia.	3.75mm	Dia.	3.75mm	Dia.	3.75mm
L	8.0mm	L	10.0mm	L	11.5mm	L	13.0mm	L	16.0mm
Ref.	91006	Ref.	91007	Ref.	91008	Ref.	91009	Ref.	91010
D1	3.85	3.85	3.85	3.85	3.85	3.85	3.85	3.85	3.85
D2	2.95	2.95	2.95	2.95	2.95	2.95	2.95	2.95	2.95
D3	3.60	3.60	3.60	3.60	3.60	3.60	3.60	3.60	3.60

PSI 4.2


Dia.	4.2mm	Dia.	4.2mm	Dia.	4.2mm	Dia.	4.2mm	Dia.	4.2mm
L	8.0mm	L	10.0mm	L	11.5mm	L	13.0mm	L	16.0mm
Ref.	91011	Ref.	91012	Ref.	91013	Ref.	91014	Ref.	91015
D1	4.20	4.20	4.20	4.20	4.20	4.20	4.20	4.20	4.20
D2	3.00	3.00	3.00	3.00	3.00	3.00	3.00	3.00	3.00
D3	3.85	3.85	3.85	3.85	3.85	3.85	3.85	3.85	3.85

PSI 5.0

Dia.	5.0mm	Dia.	5.0mm	Dia.	5.0mm	Dia.	5.0mm	Dia.	5.0mm
L	8.0mm	L	10.0mm	L	11.5mm	L	13.0mm	L	16.0mm
Ref.	91016	Ref.	91017	Ref.	91018	Ref.	91019	Ref.	91020
D1	5.00	5.00	5.00	5.00	5.00	5.00	5.00	5.00	5.00
D2	3.70	3.70	3.70	3.70	3.70	3.70	3.70	3.70	3.70
D3	3.85	3.85	3.85	3.85	3.85	3.85	3.85	3.85	3.85



■ PSI
Precise Spiral Implant



The decision on how deep to drill and what drill to use, out of those specified in the above table, is based on the accumulated experience of many dentists

Nevertheless, in cases where bone density deviates from the average, necessary adjustments should be made, relying on the experience of the dentist carrying out the procedure

Please note that in cases of soft bone, the drilling diameter is smaller than the implant's diameter, and generally the harder the bone, the larger the drilling diameter should be



PST

■ PST

Precise Spiral Top

Based on PSI's characteristics and advantages, this implant is at the top of the dental implant world today in terms of aesthetics.





■ PST Precise Spiral Top



PST – Precise Spiral Top

Topnotch Aesthetics and Stabilization


This new and advanced generation of implants is distinguished by the inverted conical design of its upper part and the differential fixation channels. This structure provides much better fixation to the bone, allows easy insertion and higher stabilization and offers greater aesthetics.

The PST implant is available in two diameters only – 4.2mm and 5mm, matching standard abutments.



Internal Hex connection

is extremely precise, being the platform for all diameters that enables a simple restoration process.



Upper Part

Surface is rough towards the top of the implant, reducing bone resorption and crestal stress for better load distribution.

2 Micro Rings



Implant Structure

Is conical, with an inverted cone in the upper part. Beyond its initial stability and easy insertion, this structure offers a clear aesthetic advantage.

Implant Threads

Are double-threaded (2x2.4mm) for efficient implantation. The threads in the upper part are thicker square, in the middle part thinner square and in the lower part V-shaped to enable fast insertion, bone condensing and self tapping.




Lower part

Has narrow, deep and sharp apical starting threads that enable dual entrance, self tapping and self drilling, as well as easy bone penetration to prevent damage to the anatomical structure.






■ PST



The PST implant is available in two diameters only - 4.2mm and 5mm - to be matched with the standard abutments.



PST 4.2


Dia.	4.2mm	Dia.	4.2mm	Dia.	4.2mm	Dia.	4.2mm	Dia.	4.2mm
L	8.0mm	L	10.0mm	L	11.5mm	L	13.0mm	L	16.0mm
Ref.	91075	Ref.	91076	Ref.	91077	Ref.	91078	Ref.	91079
D1	4.20		4.20		4.20		4.20		4.20
D2	3.00		3.00		3.00		3.00		3.00
D3	3.85		3.85		3.85		3.85		3.85
D4	3.85		3.85		3.85		3.85		3.85

PST 5.0

Dia.	5.0mm	Dia.	5.0mm	Dia.	5.0mm	Dia.	5.0mm	Dia.	5.0mm
L	8.0mm	L	10.0mm	L	11.5mm	L	13.0mm	L	16.0mm
Ref.	91080	Ref.	91081	Ref.	91082	Ref.	91083	Ref.	91084
D1	5.00		5.00		5.00		5.00		5.00
D2	3.70		3.70		3.70		3.70		3.70
D3	3.85		3.85		3.85		3.85		3.85
D4	4.40		4.40		4.40		4.40		4.40



■ PST
Precise Spiral Top



The decision on how deep to drill and what drill to use, out of those specified in the above table, is based on the accumulated experience of many dentists

Nevertheless, in cases where bone density deviates from the average, necessary adjustments should be made, relying on the experience of the dentist carrying out the procedure

Please note that in cases of soft bone, the drilling diameter is smaller than the implant's diameter, and generally the harder the bone, the larger the drilling diameter should be



PSS

■ PSS
Precise Spiral Standard

An improved iteration of the standard spiral conical implant with superior finishing, which greatly outperforms competing products.



SECTION A-A
SCALE 7:1

NOTE:
①=Critical Dimensions

NORIS Medical

IMPLANT HO



■ PSS Precise Spiral Standard



PSS – Precise Spiral Standard

The geometry and features of our PSS dental implant family, along with our meticulous manufacturing process, result in implants that ensure maximum precision, smooth insertion during the procedure and a more comfortable, natural healing process for the patient. The PSS family employs classic implant design and surface treatments to achieve a larger bone-to-implant surface area contact.



Internal Hex

Is a high-precision structure that creates the perfect implant-abutment connection and results in a simple restoration process.



Upper Part

Has rough surface through the top which reduce bone resorption and crestal stress while improving load distribution.

2 Micro Rings



Implant Structure

Is conical in nature to achieve primary stabilization and easy insertion.

Implant Threads

Are double-threaded (2x2.4mm) for efficient implantation. The threads in the upper part are thicker square, in the middle part thinner square and in the lower part V-shaped to enable fast insertion, bone condensing and self tapping.




Lower Part

Has narrow, deep and sharp apical starting threads that enable dual entrance, self tapping and self drilling, as well as easy bone penetration to prevent damage to the anatomical structure.






■ PSS



D3

D1



PSS 3.3

Dia.	3.3mm	Dia.	3.3mm	Dia.	3.3mm	Dia.	3.3mm	Dia.	3.3mm
L	8.0mm	L	10.0mm	L	11.5mm	L	13.0mm	L	16.0mm
Ref.	91050	Ref.	91051	Ref.	91052	Ref.	91053	Ref.	91054
D1	3.70	3.70	3.70	3.70	3.70	3.70	3.70	3.70	3.70
D2	2.95	2.95	2.95	2.95	2.95	2.95	2.95	2.95	2.95
D3	3.50	3.50	3.50	3.50	3.50	3.50	3.50	3.50	3.50

PSS 3.75

Dia.	3.75mm	Dia.	3.75mm	Dia.	3.75mm	Dia.	3.75mm	Dia.	3.75mm
L	8.0mm	L	10.0mm	L	11.5mm	L	13.0mm	L	16.0mm
Ref.	91055	Ref.	91056	Ref.	91057	Ref.	91058	Ref.	91059
D1	3.85	3.85	3.85	3.85	3.85	3.85	3.85	3.85	3.85
D2	2.95	2.95	2.95	2.95	2.95	2.95	2.95	2.95	2.95
D3	3.60	3.60	3.60	3.60	3.60	3.60	3.60	3.60	3.60

PSS 4.2

Dia.	4.2mm	Dia.	4.2mm	Dia.	4.2mm	Dia.	4.2mm	Dia.	4.2mm
L	8.0mm	L	10.0mm	L	11.5mm	L	13.0mm	L	16.0mm
Ref.	91060	Ref.	91061	Ref.	91062	Ref.	91063	Ref.	91064
D1	4.20	4.20	4.20	4.20	4.20	4.20	4.20	4.20	4.20
D2	3.00	3.00	3.00	3.00	3.00	3.00	3.00	3.00	3.00
D3	3.85	3.85	3.85	3.85	3.85	3.85	3.85	3.85	3.85


PSS 5.0

Dia.	5.0mm	Dia.	5.0mm	Dia.	5.0mm	Dia.	5.0mm	Dia.	5.0mm
L	8.0mm	L	10.0mm	L	11.5mm	L	13.0mm	L	16.0mm
Ref.	91065	Ref.	91066	Ref.	91067	Ref.	91068	Ref.	91069
D1	5.00	5.00	5.00	5.00	5.00	5.00	5.00	5.00	5.00
D2	3.70	3.70	3.70	3.70	3.70	3.70	3.70	3.70	3.70
D3	3.85	3.85	3.85	3.85	3.85	3.85	3.85	3.85	3.85



■ PSS

Precise Spiral Standard



The decision on how deep to drill and what drill to use, out of those specified in the above table, is based on the accumulated experience of many dentists

Nevertheless, in cases where bone density deviates from the average, necessary adjustments should be made, relying on the experience of the dentist carrying out the procedure

Please note that in cases of soft bone, the drilling diameter is smaller than the implant's diameter, and generally the harder the bone, the larger the drilling diameter should be



PHI

■ PHI

Precise Helical Implant

Our classical, mainstream, implant is an indispensable element of every implant system available in the market. The PHI promotes bone condensing and creates overall stability, it is ideal for hard bone but suitable for common cases too, and it enables small changes in the implantation direction.



SECTION A-A
SCALE 7:1

NOTE:

①=Critical Dimensions

NORIS Medical

IMPLANT HO...



■ PHI

Precise Helical Implant



The advantages of using PHI:

- Promotes bone condensing and creates overall stability
- Ideal for hard bone but can be used in common cases as well
- Enables small changes in the implantation direction

These advantages ensure easy insertion, excellent stability, self drilling and self tapping. In addition, during the implantation procedure, the PHI can be taken out slightly by turning it in the opposite direction in order to change the angle of penetration and reinsert it.



Improved Internal Hex

Has one platform for all diameters with the perfect implant-abutment connection for easy restoration process.



Upper Part

With its micro rings has a rough surface all the way to the top, which reduces crestal resorption for better load distribution and less stress on the cortical bone.

2 Micro Rings



Implant Structure

Is generally conical; while in the upper and middle parts the surface is cylindrical, in the apical part the surface tapers to achieve initial stability and easy insertion.

Implant Threads

Have steps of 2x1.2mm. The threads in the upper part are thicker square, in the middle part thinner square and in the lower part V-shaped to enable primary stabilization, easy insertion, bone condensing, self drilling and self tapping.





Lower Part

Has sharp apical threads with straight apical outlines. It is self tapping and self drilling, enabling easy penetration thereby preventing damage to the anatomical structure.





■ PHI



PHI 3.3

Dia.	3.3mm	Dia.	3.3mm	Dia.	3.3mm	Dia.	3.3mm	Dia.	3.3mm
L	8.0mm	L	10.0mm	L	11.5mm	L	13.0mm	L	16.0mm
Ref.	91025	Ref.	91026	Ref.	91027	Ref.	91028	Ref.	91029
D1	3.70	3.70	3.70	3.70	3.70	3.70	3.70	3.70	3.70
D2	2.60	2.60	2.60	2.60	2.60	2.60	2.60	2.60	2.60
D3	3.60	3.60	3.60	3.60	3.60	3.60	3.60	3.60	3.60

PHI 3.75

Dia.	3.75mm	Dia.	3.75mm	Dia.	3.75mm	Dia.	3.75mm	Dia.	3.75mm
L	8.0mm	L	10.0mm	L	11.5mm	L	13.0mm	L	16.0mm
Ref.	91030	Ref.	91031	Ref.	91032	Ref.	91033	Ref.	91034
D1	3.75	3.75	3.75	3.75	3.75	3.75	3.75	3.75	3.75
D2	3.00	3.00	3.00	3.00	3.00	3.00	3.00	3.00	3.00
D3	3.60	3.60	3.60	3.60	3.60	3.60	3.60	3.60	3.60

PHI 4.2

Dia.	4.2mm	Dia.	4.2mm	Dia.	4.2mm	Dia.	4.2mm	Dia.	4.2mm
L	8.0mm	L	10.0mm	L	11.5mm	L	13.0mm	L	16.0mm
Ref.	91035	Ref.	91036	Ref.	91037	Ref.	91038	Ref.	91039
D1	4.20	4.20	4.20	4.20	4.20	4.20	4.20	4.20	4.20
D2	3.00	3.00	3.00	3.00	3.00	3.00	3.00	3.00	3.00
D3	3.85	3.85	3.85	3.85	3.85	3.85	3.85	3.85	3.85


PHI 5.0

Dia.	5.0mm	Dia.	5.0mm	Dia.	5.0mm	Dia.	5.0mm	Dia.	5.0mm
L	8.0mm	L	10.0mm	L	11.5mm	L	13.0mm	L	16.0mm
Ref.	91040	Ref.	91041	Ref.	91042	Ref.	91043	Ref.	91044
D1	5.00	5.00	5.00	5.00	5.00	5.00	5.00	5.00	5.00
D2	4.00	4.00	4.00	4.00	4.00	4.00	4.00	4.00	4.00
D3	3.85	3.85	3.85	3.85	3.85	3.85	3.85	3.85	3.85



■ PHI

Precise Helical Implant




The decision on how deep to drill and what drill to use, out of those specified in the above table, is based on the accumulated experience of many dentists

Nevertheless, in cases where bone density deviates from the average, necessary adjustments should be made, relying on the experience of the dentist carrying out the procedure

Please note that in cases of soft bone, the drilling diameter is smaller than the implant's diameter, and generally the harder the bone, the larger the drilling diameter should be



Drills & Tools



SECTION A-A
SCALE 7:1

NOTE:
①=Critical Dimensions

NORIS Medical

IMPLANT HO



Surgical Drills

Precise's drills are available in variety of diameters. Precise's surgical drills are made of surgical stainless steel. The drills are color coded for easy diameter identification and Laser-marks on depth grooves for depth reference during surgery.

Surgical Drill
Stainless steel
Standard



○	P-DRLS 2
Ref.	94016
D	2 mm



●	P-DRLS 2.8
Ref.	94017
D	2.8 mm



●	P-DRLS 3.2
Ref.	94018
D	3.2 mm



●	P-DRLS 3.65
Ref.	94019
D	3.65 mm



●	P-DRLS 4.1
Ref.	94020
D	4.1 mm



●	P-DRLS 4.5
Ref.	94021
D	4.5 mm



●	P-DRLS 4.8
Ref.	94059
D	4.8 mm




Surgical Drill
Stainless steel
Short

		
○ P-DRLSS 2	● P-DRLSS 2.8	● P-DRLSS 3.2
Ref. 94022	Ref. 94023	Ref. 94024
D 2 mm	D 2.8 mm	D 3.2 mm

			
● P-DRLSS 3.65	● P-DRLSS 4.1	● P-DRLSS 4.5	● P-DRLSS 4.8
Ref. 94025	Ref. 94026	Ref. 94027	Ref. 94061
D 3.65 mm	D 4.1 mm	D 4.5 mm	D 4.8 mm

Other Drills



			
P-DRLE	P-MDRL 01	P-CTS 3	P-CTS 5
Ref. 94038	Ref. 4062	Ref. 94039	Ref. 94066
Extention drill	Marking drill	Countersink	Countersink

D = 3.75 x 4.2








D = 5.0 x 6.0



Surgical Accessories

Ratchet Wrenches	P-RAT	Ref.	94001	
Torque Ratchet	P-TRAT	Ref.	94002	
Screw Driver	P-SSD	Ref.	94053	

Implant Insertion Drivers

Internal Hex Insertion Drivers (stainless steel)				
	P-SRDH L	P-SRDH	P-SRDH S	
	Ref. 94003	Ref. 94028	Ref. 94029	
Internal Hex Contra-Angle Drivers (stainless steel)				
	P-IHM 1	P-IHML 2	P-IHM 3	P-IHM 4
	Ref. 94043	94044	94045	94046
	D 2.5/1.25mm	2.5mm	2.5/1.25mm	2.5mm
Internal hex motor mount	L 25mm	L 25mm	L 18mm	L 18mm





Prosthetics Insertion Drivers



	P-SRD	P-SRD L (long)		P-HSRD L (long)	P-HSRD		P-CAD 05
	Screw Driver	Screw Driver		Hand Screw Driver	Screw Driver		Contra-Angle Motor Mount
Ref.	94032	94048	Ref.	94030	94031	Ref.	94047

Parallel and Depth Guides



Depth Guide 16mm	P-DM 01		Depth Guide 10mm	P-DM 02		Parallel Guide	P-PDM
Ref.	94049		Ref.	94050		Ref.	94051

Implant Depth Prob - IDG



	P-DP 01
Ref.	94033



■ Impression Transfer



Open Tray Transfer

P-TOT 01

Ref. 92252



Open Tray Transfer (short)

P-TOT 02

Ref. 92253



Closed Tray Transfer

P-TCT 01

Ref. 92250



Closed Tray Transfer (short)

P-TCT 02

Ref. 92251



Long Screw for open Tray

P-TRS 01

Ref. 93003



Screw for closed Tray

P-TRS 02

Ref. 93004

■ Implant Analog



Analog

P-ANL 01

Ref. 92255



Analog for Wide Platform 5.0

P-ANL 02

Ref. 92256



Prosthetics & Abutments



SECTION A-A
SCALE 7:1

NOTE:
①=Critical Dimensions

REV.	DATE	DESCRIPTION	BY	CHKD.
1	09/02/11	ISSUE FOR PRODUCTION
2	09/02/11
3	09/02/11







NORIS Medical

IMPLANT HO...









■ Healing abutments

Standard Platform Healing Abutment

						
	P-HC 20	P-HC 30	P-HC 40	P-HC 50	P-HC 60	P-HC 70
Ref.	2220	2221	2222	2223	2224	2225
D	4.6mm	4.6mm	4.6mm	4.6mm	4.6mm	4.6mm
L	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm	7.0mm

Narrow Platform Healing Abutment


						
	P-HCN 20	P-HCN 30	P-HCN 40	P-HCN 50	P-HCN 60	P-HCN 70
Ref.	2226	2227	2228	2229	2230	2231
D	3.85mm	3.85mm	3.85mm	3.85mm	3.85mm	3.85mm
L	2.0mm	3.0mm	4.0mm	5.0mm	6.0mm	7.0mm



Cover Screw


P-CVS

Ref. 93002






Prosthetic Abutments



Titanium
Straight
Abutments




	P-ABTS 0002	P-ABTS 0003	P-ABTS 0004	P-ABTS 0005
Ref.	92004	92005	92006	92007
H	8.50mm	9.50mm	10.5mm	11.5mm
D	4.50mm	4.50mm	4.50mm	4.50mm
A	1.00mm	2.00mm	3.00mm	4.00mm



Titanium
Straight
Anatomic
Abutments



	P-ABTA 0012	P-ABTA 0013	P-ABTA 0014	P-ABTA 0015
Ref.	92014	92015	92016	92017
H	8.5mm	9.5mm	10.5mm	11.5mm
D	4.90mm	4.90mm	4.90mm	4.90mm
A	0.00mm	1.00mm	2.00mm	3.00mm
B	1.00mm	2.00mm	3.00mm	4.00mm




Titanium
Straight
Abutments



	P-ABT 0008	P-ABT 0011	Abutment Screw	P-COS
Ref.	92030	92031		Ref. 93001
H	8.50mm	11.5mm		
D	4.50mm	4.50mm		
A	0.80mm	0.80mm		




Prosthetic Abutments



Titanium
15° Angular
Abutments




	P-ABT 1502	P-ABT 1503	P-ABT 1504
Ref.	92011	92012	92013
H	7.50mm	9.50mm	11.5mm
D	4.50mm	4.50mm	4.50mm
A	1.30mm	1.30mm	1.30mm



Titanium
15° Angular
Anatomic
Abutments



	P-ABTA 1511	P-ABTA 1512	P-ABTA 1513	P-ABTA 1514
Ref.	92021	92022	92023	92029
H	9.50mm	10.5mm	11.5mm	12.5mm
D	5.10mm	5.10mm	5.10mm	5.10mm
A	1.00mm	2.00mm	3.00mm	4.00mm
B	2.00mm	3.00mm	4.00mm	5.00mm




Titanium
25° Angular
Abutments



	P-ABT 2502	P-ABT 2503	P-ABT 2504
Ref.	92008	92009	92010
H	7.00mm	9.50mm	11.5mm
D	4.50mm	4.50mm	4.50mm
A	1.3mm	1.3mm	1.3mm
B	0.30mm	0.30mm	0.30mm



Prosthetic Abutments



Titanium
25° Angular
Anatomic
Abutments



	P-ABTA 2511	P-ABTA 2512	P-ABTA 2513	P-ABTA 2514
Ref.	92018	92019	92020	92028
H	9.50mm	10.5mm	11.5mm	12.5mm
D	5.10mm	5.10mm	5.10mm	5.10mm
A	1.00mm	2.00mm	3.00mm	4.00mm
B	2.00mm	3.00mm	4.00mm	5.00mm

Prosthetics-Overdenture

Straight
Ball
Attachments



	P-BA 0.5	P-BA 2	P-BA 3	P-BA 4	P-BA 5	P-BA 6
Ref.	2258	2259	2260	2261	2262	2263
D	0.5	2.0	3.0	4.0	5.0	6.0

Nylon Caps & Metal House



Stainless
Steel Housing

	P-HSB
Ref.	92236



Nylon Cap

	P-CFB 01
Ref.	92234



Nylon Cap
With Titanium

	P-CFBT 02
Ref.	92235

Plastic Abutments



Material
plastic

	P-PAB 01	P-PAB 02
Ref.	2232	2233



Precise Implant Systems Ltd.
P.O.B 45469
Tel Aviv
Postcode 6145401