

HemoSplit[®]
LONG-TERM HEMODIALYSIS CATHETER

HemoSplit[®] XK
LONG-TERM HEMODIALYSIS CATHETER

With
Air Guard[™]
VALVED INTRODUCER



Extreme³ Access[™]

BARD

ACCESS SYSTEMS

EXCEPTIONAL PERFORMANCE

EASY INSERTION

EFFECTIVE TECHNOLOGY

**Equipment You Can Count on...
for Extreme Performance in Long-Term Dialysis and
Improved Protection from Blood Loss and Air Embolism.**

**Extreme
Access™**

HemoSplit®
LONG-TERM HEMODIALYSIS CATHETER

Exceptional Performance

- **Carbothane catheter material** – affords strength for longevity and softness for flexibility and patient comfort.
- **360° multiple sideholes** – reduce the risk of catheter occlusion by the vessel wall, thus reducing the risk of arterial insufficiency.
- **Large lumens and non-restrictive tip design** – enable flow rates as high as 500ml/min.

Easy Insertion

- **Exceptional kink resistance** – simplifies insertion by allowing greater flexibility in tunnel location.
- **Guidewire passage/hole** – on the venous tip facilitates threading for easy over-the-wire insertion.
- **Malleable tunneler, PTFE introducer and stiff guidewire** – facilitate catheter insertion.

**AirGuard™ Valved Introducer Improved
Protection from Air Embolism and Blood Loss**

- **0.15cc/second air leakage*** – Most effective valved introducer for preventing air embolism.
- **Reduction in blood loss** – A safer procedure when compared to non-valved introducers.



Integrated Valve

* At -12mm Hg Vacuum Pressure

Please consult product labels and inserts for any indications, contraindications, hazards, warnings, cautions and instructions for use.

Bard Access Systems, Inc.

United States


Bard Access Systems
5425 West Amelia Earhart Drive
Salt Lake City, Utah 84116
801-595-0700

Customer Service: 800-545-0890
Clinical Information: 800-443-3385
www.bardaccess.com

Canada

Bard Canada
2345 Stanfield Road
Mississauga, Ontario L4Y3Y3

Customer Service:
Western Canada 800-268-2862
Eastern Canada 800-387-7851
Ontario 800-387-9473



HemoSplit^{XX}
LONG-TERM HEMODIALYSIS CATHETER

**Optimized Performance for
Difficult Patients**

HemoSplit XK catheter offers the same advantages as the **HemoSplit** catheter with the following features:

- **16 Fr. catheter**—allows for flow rates of 500mL/min with lower pressure than the 14.5 Fr. **HemoSplit** Catheter.
- **Increased longevity**—larger inner-lumen area allows for flows of 500 mL/min *even when occluded 15%*.

14.5 Fr. **HemoSplit**[®]
Carbothane Shaft

16 Fr. **HemoSplit^{XX}**
Carbothane Shaft

Separated
Split Tips

HemoSplit XK 16 Fr. Product Codes		HemoSplit 14.5 Fr. Product Codes			
Standard Kit	Microintroducer Kit	Standard Kit	Microintroducer Kit	Tip to Cuff	Tip to Bifurcation
Straight					
5683690	5693690	5733150		15cm	20cm
5683730	5693730	5733690	5743690	19cm	24cm
5683270	5693270	5733730	5743730	23cm	28cm
		5733270	5743270	27cm	32cm
5683310	5693310	5733310	5743310	31cm	36cm
5683350		5733350	5743350	35cm	40cm
5684420		5734420	5744420	42cm	47cm